

Labycare-Edridge Green Lantern (SKSS-MD2044)



Edridge Green Lantern Manufacturer by LABYCARE

Description :

Edridge Green lantern is a diagnostic device used to test color vision, especially for detecting red-green color blindness. It simulates signal lights viewed at a distance. Commonly used in ophthalmology and occupational screening, it helps assess an individual's ability to correctly identify colored signals.

Specifications :

Type: Color vision testing lantern

Light Source: Red, green, white filters

Operation: Manual/electric illumination

Viewing Distance: Adjustable

Use: Color blindness and signal recognition test

Application: Clinical and occupational screening

labycare

www.labycare.com | Email: labycare@gmail.com