


Labycare-Snellen"s Charts (SKSS-MD2045)



Snellen"s Charts Manufacturer by LABYCARE

Description : Snellen's chart is a standard eye testing chart used to measure visual acuity. It consists of rows of letters decreasing in size, viewed from a fixed distance. Widely used in clinics and hospitals, it helps assess clarity of vision and detect refractive errors accurately.

Specifications :

Type: Visual acuity chart

Material: Paper / plastic laminated

Design: Letters of decreasing size

Standard Distance: 6 meters (20 feet)

Illumination: Backlit / wall-mounted

Use: Vision testing and screening

labycare

www.labycare.com | Email: labycare@gmail.com