

Labycare-Maddox Rod (SKSS-MD2046)



Maddox Rod Manufacturer by LABYCARE

Description :

Maddox rod is an ophthalmic instrument used to detect ocular misalignment and diagnose strabismus or diplopia. It consists of a series of parallel cylindrical lenses that convert a point of light into a line. It is commonly used in eye examinations for assessing binocular vision and alignment.

Specifications :

Type: Ophthalmic diagnostic lens
Material: Glass/plastic cylindrical rods
Color: Red (commonly)
Design: Parallel rod structure
Use: Ocular alignment and diplopia testing

Application: Clinical eye examination

labycare

www.labycare.com | Email: labycare@gmail.com