

Labycare-Analgesiometer (SKSS-MD2049)



Analgesiometer Manufacturer by LABYCARE

Description :

Analgesiometer is a laboratory instrument used to measure pain threshold and analgesic effects in experimental subjects. It applies controlled thermal or mechanical stimulus to evaluate response. Widely used in pharmacology research, it helps assess effectiveness of pain-relieving drugs and study nociceptive mechanisms accurately.

Specifications :

Type: Pain threshold measurement device
Method: Thermal (hot plate/tail flick) or mechanical
Temperature Range: Adjustable (for thermal type)
Timer: Digital/analog response timer
Controls: Variable intensity settings
Display: Digital/analog readout
Power Supply: 220–230V AC, 50/60 Hz

labycare

www.labycare.com | Email: labycare@gmail.com