

## Labycare-Digital Photo Actometer (SKSS-MD2052)



Digital Photo Actometer Manufacturer by LABYCARE

### **Description :**

Digital photo actometer is a laboratory instrument used to measure locomotor activity of experimental animals. It uses infrared beams to detect movement and record activity counts. Widely used in pharmacology research, it helps evaluate stimulant or depressant effects of drugs on central nervous system activity.

### **Specifications :**

Type: Locomotor activity measurement device

Detection: Infrared beam sensors

Display: Digital activity counter

Chamber: Enclosed activity cage

Sensitivity: Adjustable

Timer: Built-in timer function

Power Supply: 220–230V AC

Use: CNS and behavioral studies

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)