

## Labycare-Rotarod Apparatus Various (Product Cat # SKSS-MD2054)



Rotarod Apparatus Various Manufacturer by LABYCARE

### **Description :**

Rotarod apparatus is a laboratory device used to evaluate motor coordination, balance, and grip strength in experimental animals. It measures the time an animal remains on a rotating rod. Commonly used in pharmacology studies, it helps assess effects of drugs on motor function and neurological performance.

### **Specifications :**

Type: Motor coordination testing device  
Rod: Rotating horizontal rod  
Speed: Adjustable (constant/accelerating)  
Lanes: Single or multiple lanes  
Timer: Fall-off time recording  
Display: Digital/analog  
Power Supply: 220–230V AC

Use: Neurological and pharmacological studies

---

labycare

www.labycare.com | Email: labycare@gmail.com