

Labycare-Physiograph Various (SKSS-MD2055)



Physiograph Various Manufacturer by LABYCARE

Description :

Physiograph is a multi-channel recording device used to measure and record physiological signals such as ECG, respiration, muscle activity, and blood pressure. It converts biological signals into graphical data for analysis. Widely used in laboratories, it supports accurate monitoring and study of various body functions.

Specifications :

Type: Multi-channel physiological recorder
Channels: 2–8 or more (model dependent)
Parameters: ECG, respiration, EMG, BP
Display: Digital/graphical output
Input: Sensor/electrode compatible
Data: Real-time recording and storage
Power Supply: 220–230V AC

Use: Physiology and research applications

labycare

www.labycare.com | Email: labycare@gmail.com