

## Labycare-Clamps Various (SKSS-MD2056)



Clamps Various Manufacturer by LABYCARE

### **Description :**

Clamps are laboratory accessories used to hold and secure apparatus such as burettes, flasks, and tubes during experiments. They provide stability and proper positioning of equipment. Commonly used with stands, clamps ensure safe handling and accurate experimental setup in laboratories and medical applications.

### **Specifications :**

Type: Burette clamp / boss head / 3-prong clamp

Material: Mild steel / stainless steel

Finish: Chrome plated / powder coated

Adjustment: Screw tightening mechanism

Compatibility: Retort stand mounting

Use: Holding laboratory apparatus

---

labycare

www.labycare.com | Email: labycare@gmail.com