

Labycare-Travelling Microscope (SKSS-MD2064)



Travelling Microscope Manufacturer by LABYCARE

Description :

The Travelling Microscope is a precision instrument used to measure very small distances accurately. It consists of a microscope mounted on a movable scale with vernier readings. It is widely used in laboratories for optical experiments, refractive index determination, and precise dimensional measurements of small objects.

Specifications :

Base: Heavy cast iron/metal base with leveling screws

Movement: Smooth horizontal and vertical sliding arrangement

Travel Range: 150–220 mm (horizontal), 140–150 mm (vertical)

Least Count: 0.01 mm or 0.02 mm

Scale: Engraved metal scale with vernier

Eyepiece: 10× with crosshair

Objective Lens: Achromatic lens

labycare

www.labycare.com | Email: labycare@gmail.com