

Labycare-Bunsen Burner with Stopcock For Labware Biology Chemistry (SKSS-MD2069)



Bunsen Burner with Stopcock For Labware Biology Chemistry

Description :

The Bunsen Burner is a common laboratory instrument used for heating, sterilization, and combustion. Equipped with a stopcock, it allows precise control of gas flow and flame intensity. Widely used in biology and chemistry labs, it ensures safe, efficient, and adjustable heating for various experimental procedures.

Specifications :

- Adjustable stopcock for precise gas control
- Stable flame with air inlet regulation
- Durable metal construction
- Suitable for multiple gas types

Easy to operate and maintain

labycare

www.labycare.com | Email: labycare@gmail.com