

Labycare-Spirit Lamp (SKS-SE01)



Spirit Lamp

Description :

Spirit lamp is a laboratory heating device that uses alcohol as fuel to produce a clean flame. It is commonly used for gentle heating, sterilization, and flame tests. Compact and easy to use, it is suitable for educational and laboratory applications requiring low-temperature heating.

Specifications :

Type: Alcohol-based burner

Material: Glass / metal body

Fuel: Ethanol / methanol

Capacity: Approx. 100–250 ml

Wick: Cotton wick with holder

Flame: Adjustable via wick length

Cap: Metal cap for extinguishing