

Labycare-Retort Stand (SKS-SE02)



Retort Stand

Description :

Retort stand is a laboratory support apparatus used to hold and position equipment such as clamps, burettes, and flasks. It consists of a heavy base and vertical rod, providing stability. Widely used in laboratories, it ensures safe and organized setup during experiments and procedures.

Specifications :

Material: Cast iron base with steel rod

Base: stable and non-slip

Rod Height: Approx. 45–100 cm

Compatibility: Clamps, boss head attachments

Finish: Powder coated / chrome plated

Use: Supporting laboratory apparatus

labycare

www.labycare.com | Email: labycare@gmail.com