

Labycare-Test Tube Holder (SKSS-SE08)



Test Tube Holder

Description :

A test tube holder is a handheld laboratory tool designed to safely grip and hold test tubes during heating or chemical reactions. It protects the user from burns and chemical exposure. Typically made of wood or metal, it provides a firm grip and ease of handling in laboratory experiments.

Specifications :

Material: Wood / Stainless Steel / Mild Steel

Length: 150 mm to 200 mm (approx.)

Grip Type: Spring-loaded or clamp-style

Handle: Insulated (in wooden models)

Corrosion Resistance: High (for metal types)

Heat Resistance: Suitable for high-temperature handling

labycare

www.labycare.com | Email: labycare@gmail.com