

Labycare-Vernier Caliper (Digital) (SKSS-SE030)



Vernier Caliper (Digital)

Description :

Digital Vernier Caliper is a precision measuring instrument used for accurate measurement of internal, external, depth, and step dimensions. It has an LCD display for easy reading, zero setting, and unit conversion. Made from stainless steel, it ensures durability, high accuracy, and reliable performance in laboratory and industrial applications.

Specifications :

Measuring Range: 0–150 mm (6 inch)

Resolution: 0.01 mm / 0.0005 inch

Accuracy: ± 0.02 mm

Material: Hardened stainless steel

Display: Digital LCD

Power Supply: Button cell battery

Functions: Internal, external, depth, and step measurement

labycare

www.labycare.com | Email: labycare@gmail.com