

Labycare-Micrometer/Screw Gauge (Digital) (SKSS-SE031)



Description :

Digital Micrometer (Screw Gauge) is a high-precision instrument used to measure small dimensions like thickness and diameter accurately. It features an LCD display for clear readings, ratchet stop for uniform pressure, and zero setting. It ensures reliable, precise measurements and is widely used in laboratories, engineering, and quality control work.

Specifications :

Measuring Range: 0–25 mm (standard)

Resolution: 0.001 mm / 0.00005 inch

Accuracy: ± 0.002 mm

Display: Digital LCD screen

Material: Hardened steel frame with carbide measuring faces

Power Supply: Button cell battery

Functions: External measurement, zero setting, mm/inch conversion

labycare

www.labycare.com | Email: labycare@gmail.com