

Labycare-Analytical Weight Box (SKSS-SE033)



Analytical Weight Box

Description :

Analytical Weight Box is a precision set of standard weights used for accurate mass measurement in laboratories. It contains carefully calibrated weights of various denominations, stored in a protective case. Used with balances, it ensures reliable results, high accuracy, and is essential for calibration, testing, and scientific experiments.

Specifications :

Weight Range: 1 mg to 100 g (typical set)

Material: Stainless steel / brass (polished)

Accuracy Class: E2 / F1 (high precision)

Type: Knob type and sheet type weights

Box: Wooden or plastic case with cushioned lining

Accessories: Forceps for handling small weights

Features: Corrosion resistant, precisely calibrated, durable

labycare

www.labycare.com | Email: labycare@gmail.com