

Labycare-SLOTTED WEIGHT BRASS/IRON (SKSS-SE034)



Description :

Slotted Weights are precision masses made of brass or iron with a central slot for easy placement on a hanger. Used in physics experiments, they help study force, mass, and motion. These weights are durable, accurately calibrated, and commonly used in laboratories for mechanical and educational applications.

Specifications :

Material: Brass / Cast iron

Available Weights: 10 g, 20 g, 50 g, 100 g, 200 g (or set)

Design: Slotted with central opening for hanger

Accuracy: Standard laboratory precision

Finish: Polished (brass) / painted (iron)

Compatibility: Suitable for weight hanger/stand

Application: Physics experiments, force and motion studies