

Labycare-UNIVERSAL LENS HOLDER (SKSS-SE037)



UNIVERSAL LENS HOLDER

Description :

Universal Lens Holder is a laboratory apparatus used to hold optical lenses securely during experiments. It is designed to fit lenses of various sizes and is mounted on a stand for stability. Commonly used in optics experiments, it ensures proper alignment and accurate positioning for studying light and image formation.

Specifications :

Material: Metal frame with durable finish
Design: Adjustable holder for different lens diameters
Mounting: Compatible with optical bench or stand
Clamp Type: Spring or screw-based grip

Height Adjustment: Provided with rod and base

labycare

www.labycare.com | Email: labycare@gmail.com