

Labycare-CHEMICAL BALANCE (SKSS-SE040)



CHEMICAL BALANCE

Description :

Chemical Balance is a precision instrument used to measure mass accurately in laboratories. It operates on a beam or digital system to compare weights. Commonly used in chemistry experiments, it ensures high accuracy and reliability for weighing chemicals, reagents, and samples in analytical and research applications.

Specifications :

Type: Mechanical (beam) / Digital balance

Capacity: 200 g to 500 g (typical)

Least Count: 0.001 g (digital) / 0.01 g (mechanical)

Accuracy: High precision measurement

Material: Metal body with glass draft shield

Display: LCD (digital type)

Calibration: Manual or automatic

labycare

www.labycare.com | Email: labycare@gmail.com