

Labycare-TICKER TAPE TIMER (SKSS-SE044)



Description :

Ticker Tape Timer is a laboratory device used to study motion and time intervals in physics experiments. It marks dots on a paper tape at regular intervals using an electric vibrator. By analyzing the spacing of dots, velocity and acceleration of moving objects can be calculated accurately.

Specifications :

Type: Electromagnetic timing device

Frequency: 50 Hz (marks 50 dots per second)

Power Supply: AC mains or low voltage supply

Output: Dotted paper tape

Material: Metal/plastic body with internal mechanism

Operation: Electric vibrator creates uniform time intervals

Application: Study of velocity, acceleration, motion analysis

labycare

www.labycare.com | Email: labycare@gmail.com