

## Labycare- TRIPLE BEAM BALANCE (SKSS-SE048)



### TRIPLE BEAM BALANCE

#### **Description :**

Triple Beam Balance is a mechanical instrument used to measure mass accurately in laboratories. It consists of three beams with sliding riders for precise measurement. Commonly used in schools and labs, it provides reliable readings without electricity and is ideal for weighing solids and small objects with good accuracy.

#### **Specifications :**

Type: Mechanical beam balance

Capacity: 610 g (with additional weights up to 2610 g)

Least Count: 0.1 g

Beams: Three graduated beams with riders

Material: Metal base with stainless steel pan

Calibration: Manual zero adjustment

Display: Mechanical scale reading

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)