

Labycare-Hospital General Bed (SKSS809M)



Hospital General Bed Manufacturer by LABYCARE

Description :

A **Labycare Hospital General Bed**, is specifically designed for non-critical hospitalized patients or others who require medical attention. These plain hospital beds have special features for the patient's comfort and well-being, as well as the convenience of health care workers.

Specifications :

- Mild steel tubular Head and Foot Bows of unequal height with vertical tube support
- Framework made up of M.S. tubes
- Top made up of oval perforated CRCA M.S. sheet
- Legs fitted with PVC stumps
- Four location for I.V. Rod
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 1900L x 910W x 600H mm

labycare

www.labycare.com | Email: labycare@gmail.com