

Labycare-Hospital Couch Bed (SKSS809F)



Hospital Couch Bed manufacturer by **Labycare** in India.

Description :

Facilitated by technically advanced process, we are involved in offering a vast range of **Labycare Hospital Couch Bed**. The offered product is praised for its attractive look and premium quality. A series of quality tests is performed on our product at the time of manufacturing by our quality controllers to ensure its quality. Moreover, this qualitative product is widely demanded in domestic as well as international markets.

Specifications :

- Overall size: 1800L x 600W x 460H mm
- Framework made up of Rectangular / Square M.S. tubes
- Also available in 5 cm from upholstery top
- Legs fitted with PVC stumps
- Finish: Pre-treated & Epoxy Powder Coated

Optional

- 5cm foam padded fixed mattress covered with rexine
- Plain Mattress

labycare

www.labycare.com | Email: labycare@gmail.com