

Labycare-Examination Bed (SKSS809N)



Examination Bed manufacturer by **Labycare** in India

Description :

Labycare Examination Bed is used to support patients during medical examinations. During these exams, doctors in offices, clinics and hospitals use an adjusting mechanism to manipulate and position the table to allow patient support, closer examination of a portion or the entire patient, and the ability to move the patient on and off the table safely.

Specifications :

Bed type: Examination bed / examination couch

Frame material: Mild steel tubular construction

Finish: Epoxy powder coated

Bed top: Foam padded with washable, waterproof upholstery

Backrest: Adjustable (manual ratchet)

Dimensions: Approx. 1800–1900 mm (L) × 600–650 mm (W) × 750–800 mm (H)

Load capacity: Approx. 150–180 kg

Stability: Anti-skid rubber feet

Maintenance: Easy to clean

labycare

www.labycare.com | Email: labycare@gmail.com