

Labycare-Semi Bed (SKSS8090)



Hospital Semi Beds manufacturer by **Labycare** in India.

Description :

Labycare Hospital Semi Beds are designed, four section perforated beds utilized for knee rest and backrest functions. The beds, operated using cranks or motors, utilize either simultaneous, central or individual braking system.

Specifications :

- ABS moulded Head and Foot Panel with locking and corner buffer
- Frame work made up of M.S. tubes
- 2 section top made up of perforated CRCA M.S. sheet
- Legs fitted with PVC stumps
- Four location for I.V. Rod
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 900W x 600H mm
- Positions obtained by Smooth Crank System with foldable handle:
- Backrest elevation

Optional

- 2 Section Mattress
 - S.S. Telescopic I.V. Rod
 - S.S. Collapsible Railing
 - Advance linear actuator system (01 Pc Motor, 01 Pc control box, 01 Handset)
 - M.S. or S.S. Head & Foot Panels
-

