

Labycare-Electric cadaver lifting tray (SKSSECLT246)



We are engaged in manufacturing and supplying of **Electric Tray Lifting** at a large scale. Our organization uses contemporary technologies and optimum range of raw materials for the designing and developing of our products.

Description :

We are engaged in manufacturing and supplying of **Electric Tray Lifting** at a large scale. Our organization uses contemporary technologies and optimum range of raw materials for the designing and developing of our products.

Specifications :

- Electric lifting trolley with double cross bar for sliding body trays with casters.
- Structure with lined articulations - Stowage of the body tray while transporting
- Electronic silent jack operates lifting and descending of the plate : 24VCC, 1x2,9Ah (about 40 lifts),
- battery supplied. Battery charger can be used from +5°C to a +40°C. Works with 230V 50/60 Hz
- Button action lifting and descending of the trolley. - Body tray holding system provided.
- Emergency stop button,"punch-stop".
- 4 rotating wheels, 100 mm diameter, 2 of them with brakes. Ball-bearing system. Can be used from
- 20°C to +60°C. Finishe : Hammer-wrought textured beige steel, water resistant.

labycare

www.labycare.com | Email: labycare@gmail.com