

Labycare-Laboratory biosafety cabinet class 2 (SKSS 395)



Description :

The Claire pure provides the highest quality technology at an attractive price. LED illumination in the work area, high-performance EC fans and the touch display with intuitive menu navigation form the technical basis for the highest functionality. This is complemented by a targeted range of options for all major applications.

Outstanding energy efficiency and low procurement and operating costs make Claire pure the ultimate in efficiency.

Specifications :

- vertical laminar air flow
- II safety class
- 2 HEPA filters with efficiency of 99.995% for MPPS
- microprocessor control system with touchscreen in Polish/English or any other language
- automatic filter compensation ensures constant and safe speed of air flow in the cabinet
- hour timer of the cabinet and UV lamp operation
- adjusted air flow
- clear and easy-to-use control panel – optical and acoustic signals
- silent-running fans,
- front window made of tempered double glass without frame, inclined at an angle of 8°, electrically raised and lowered, with a possibility to shut it completely – closing the workspace
- split worktop, made of stainless steel
- internal workspace made of stainless steel
- permanently fixed UV lamp in the upper part of the workspace, protected from accidental switching on during operation
- white fluorescent lightning
- workspace sides with window openings
- removable armrest
- frame with lockable castors
- power supply 230V/60 Hz

- *digital time counter of cabinet and inner oversee system of errors*

labycare

www.labycare.com | Email: labycare@gmail.com