

Labycare-incubator shaker (SKSSIS336)



This shaking incubator is demanded cell culturing, bacterial culturing, cell aeration, and metabolism & solubility studies. Also, offered shaking incubator is tested by our quality controllers against diverse parameters so as to dispatch a defect-free range.

Description :

We are counted among the leading manufacturer, exporter and supplier of quality assured **Shaking Incubators**. Our provided shaking incubators are appreciated by our clients for their compact design, easy installation, strength and longer service life. This shaking incubator is demanded cell culturing, bacterial culturing, cell aeration, and metabolism & solubility studies. Also, offered shaking incubator is tested by our quality controllers against diverse parameters so as to dispatch a defect-free range.

Specifications : Number of Shelves	3 Shelves
Model	SKSSIS336D
Chamber Volume	250 Litre
Temperature Range	5 to 105 degree celsius
Brand	SKSS
Voltage	230v
Frequency	50 hz
Material	ss
Automation Grade	full automation
Power	1 hp
Minimum Order Quantity	1 Unit

labycare

www.labycare.com | Email: labycare@gmail.com