

Labycare-Stability chamber (SKSSS401)



Stability Chambers adopts original imported high quality components and manufacturing process, with stable and reliable performance and three functions of temperature, humidity and illumination. It is equipped with visible light sensor, which can record and store temperature, humidity and visible light values.

Description :

labycare series comprehensive medicine Stability Chambers adopts original imported high quality components and manufacturing process, with stable and reliable performance and three functions of temperature, humidity and illumination. It is equipped with visible light sensor, which can record and store temperature, humidity and visible light values.

<u>Specifications :</u>	TEMP.RANGE (?)	HUMI.RANGE (RH)	CAPACITY (L)	INTERNAL DIMEN. (mm) W×D×H	EXTERNAL DIMEN. (mm) W×D×H
SKSSSC365	15?65	40?95%	150	600×405×620	760×830×1480
SKSSSC365A	15?65	40?95%	250	600×500×830	760×870×1650
SKSSSC365B	15?65	40?95%	400	680×600×930	840×970×1750
SKSSSC365C	15?65	40?95%	500	680×680×1080	840×1050×1850
SKSSSC365D	15 ? 65	40? 95%	800	1200×510×1300	1350×870×1880
SKSSSC365E	15 ? 65	40? 95%	1000	1400×510×1400	1550×870×1980

labycare

www.labycare.com | Email: labycare@gmail.com