

Labycare-12 hole water bath (SKSS344D)



Being a prominent Manufacturer, Supplier and Exporter in the industry, we present a wide range of **Laboratory 12 HOLE Water Bath** to our customers. The quality we offer is unmatched. We take special care of quality by matching it against the highest parameters of quality set by the industry. We provide it in safe packaging at affordable rates. The range of offer enjoys a long life and lustrous design

Description :

- This water baths are ideal for general purpose application. Bath fluids can be controlled at temperatures as high as 100oC (60oC without cover) with ± 0.20 precision and uniformity. Bath temperature can be read from mercury thermometer (included) with a range from -20 to 110oC. Cover and ring lids are included.
- Digital bath accommodates vessel sizes up to 2000ml and features a digital temperature controller cum timer and low water level indicator with cutoff system.

Specifications :

- -Thermostatically Water Bath is a double walled
- -The gap between the inner and outer chamber fines with special grade glass-wool insulation to avoid the heat loss.
- -Concentric rigs can be removed to accommodate containers of various sizes.
- -Temperature 5 degree Celsius above ambient to 95 degree Celsius.
- -Temperature is controlled by thermostat.

labycare

www.labycare.com | Email: labycare@gmail.com