

Labycare-Single wall water bath (SKWB345E)



Water Bath Single Wall Stainless Steel 300mm*250mm*100mm 6 holes

Description :

A rectangular water bath is a system designed for the general boiling purpose in laboratories. There are a lot of features offered by the system to make working easy in the laboratory. The body of the system is made out of thick stainless steel, top made out of chrome-plated sheet concentric rings to cover hole are provided. As the chamber is made up of stainless steel it is easy to clean. There is also a drain plug fitted for quick water drainage. The water bath system has an auto-cut heating elements are fitted to protect the same when the bath gets dry.

Specifications :

- 300mm*250mm*100mm (6 holes)(12 holes)
- Single-walled construction
- Temperature controller
- Removable concentric rings
- Stainless steel inner chamber
- Drain plug for easy cleaning
- Supplied with or without stirring facility
- Works on 220/230 volts A.C. supply

labycare

www.labycare.com | Email: labycare@gmail.com