

Labycare-Water Bath Incubator Shaker (SKSS348)



Orbital shaking motion Included spring platform accepts tubes, bottles and flask Compact design, conserves bench space Digital control and display of time, temperature, and speed.

Description :

This compact shaking water bath incubator utilizes the well proven, belt-less motor drive system from our popular -Shaker, with the addition of precision temperature control up to 80°C.

The seamless, stainless steel chamber of the water bath incubator has a 12 liter capacity and includes a universal spring-loaded platform for securing 15 & 50ml test tubes, as well as Erlenmeyer flasks and media bottles up to 1 liter. A hinged and removable lid is also included to reduce energy waste, prevent splashing and divert condensation away from samples.

<u>Specifications :</u>		
Temp. range :		Room temperature +5°C to 80°C
Capacity:		12 liters
Temp. accuracy:		±0.3°C
Temp. uniformity:		±0.2°C
Temp. increment:		0.1°C
Shaking speed ranges:		30 to 200 rpm (increments of 5 rpm)
Timer:		1 min. to 48 hr. or continuous
Dimensions : (w x d x h):		12.2 x 13.5 x 12 in./31.7 x 34.3 x 30.5 cm
Weight:		31 lbs./14kg
Electrical:		115 or 230V , 50-60Hz/300W

labycare

www.labycare.com | Email: labycare@gmail.com