

Labycare-Constant Temperature Water Bath. (SKSSWB347H)



Break through the defects of lack of water and burn prevention in existing domestic water tanks, stamping and forming stainless steel tanks at a time, with an independent water-proof and anti-air-burning function and an independent temperature-limiting safety system to ensure safe operation without accidents.

Description :

Liquid crystal display, multiple groups of data displayed on one screen, menu-type operation interface, with timing function, easy to observe and operate.

(Optional multi-stage programmable liquid crystal temperature controller)

One-time stamping forming stainless steel liner, anti-corrosion, easy to clean.

Stainless steel tip cover prevents moisture evaporation and dew from entering the test sample causing contamination.

With water shortage and anti-air heating function, the water shortage phenomenon immediately cut off the heater, and the sound and light alarm prompts the operator to ensure the safe operation of the instrument without accidents.

It can connect micro printer, 485 communication interface or USB data interface. (Optional)

Specifications :

Model Specifications	SKSSWB347H	SKSSWB347I	SKSSWB347J	SKSSWB347K
voltage	AC220V 50HZ			
input power	800W		1000W	
Temperature	RT+5 ~ 100 °C	RT+5 ~ 80 °C	RT+5 ~ 100 °C	RT+5 ~ 80 °C

control range				
Constant temperature fluctuation	±0.2 °C			
Resolution	0.1 °C			
Effective volume	11L		20L	
Liner size	300×240×200		500×300×200	
Dimensions	480×300×480		680×360×390	
Timing range	0 to 9 999min			

labycare

www.labycare.com | Email: labycare@gmail.com