

Labycare-Table Top Sieve Shaker. (SKSSTSS359C)



Table Top Sieve Shaker assists in producing uniform motion for shaking. The vibration is controlled by speed regulator. It is used for research and development, quality control of raw materials, as well as in production monitoring.

Description :

Sieve Shaker Table top type To make the process of sieving easier and quicker, electrically operated mechanical table top type sieve shakers are offered for dry sieving. These Sieve Shakers are popularly used not only in soil laboratories but also dry sieving the geotextile testing.

<u>Specifications :</u>	Model	SKSSTSS359C	SKSSTSS359D	SKSSTSS359E	SKSSTSS359F
Sieving Capacity		(10) 8" sieves (12) 5" sieves (16) 3" full-height sieves, (18) 8" half-height sieves.	(4) 12" full-height sieves (10) 12" half-height sieves	(6) 8" full-height sieves (13) 8" half-height sieves	(10) 8" and (6) 12" full-height sieves
Construction	Cast Iron, Painted Mild steel				
Motor	¼ HP				
Timer	0 - 60 Minutes				
Timer Type	Digital / Analog				
Standard	ASTM C136				
Power Supply	220 Volts				
Certifications	CE and ISO				