

Labycare-SEED GERMINATOR DOUBLE CHAMBER (SKSS-SLE2027)



- Seed Germinator Double Chamber Manufacturer, Seed Germinator Double Chamber Exporter, Heat And Refrigeration System, Seed Germinator Double Chamber, Laboratory Equipments, Germinator Double Chamber Exporter, Germinator Double Chamber Suppliers, Germinator Double Chamber Manufacturer in India.

Description :

Germinators are useful for seed testing, biological studies, forestry research works etc. With a solid see thru double walled door and a full view inner plexi-glass door enables inspection and monitoring of inner chamber specimens without disturbing the process temperature. Excellent quality magnetic gum packing door gasket for external door. User- oriented design of shelves makes you adjust each space of shelves without difficulty. Forced air circulation at triple walled back by durable coaxial blower, maintains optimum temperature uniformity and homogeneity

<u>Specifications :</u>	No. of chambers	Two (each having heating & cooling)
	Chamber size	55 x 55 x 95 cm (L x B x H) (Approx).
	Capacity	10 cuft. (Approx.) (each chamber).
	Inside Chamber	Stainless steel grade 304.
	Outside Chamber	Painted steel of 20 SWG.
	Temperature	Microprocessor based control 5°-45°C ±1°C.
	Evaporator	Finned tube
	Compressor	Hermetically sealed type of standard.
	Light	700 -1250 lux intensity of light.
	Humidity	Humidity Control Range Ambient 90%RH with Steam humidifier.
	No of trays	3: 3 adjustable SS wiremesh trays.(each chamber) Tray size 530mm x 530mm.

Temperature and humidity sensor	Accurate measurement of temperature and humidity (Swiss made rotronic hygroclip as per customer requirement).
Display	2×16 characters LCD display.
Door panel	For fitting of all electronic circuits and control panel so that maintenance is easy and control panel is at a convenient.
CE safety test standards	Meets CE safety test standard.
Heater dry and wet	CE approved.
PUF insulation	For energy saving.

labycare

www.labycare.com | Email: labycare@gmail.com