

Labycare-Examination Dental Chair (SKSS876)



The Examination Dental Chair manufacturer by **Labycare** in India.

Description :

Labycare Examination Dental Chair is a specially designed medical device intended to support a patient's body when a dental procedure is being performed on them. _____

Equipped with advanced accessories such as an intraoral camera, LED curing unit, and digital X-ray system, it supports accurate diagnosis and efficient treatment, making it ideal for modern dental clinics and hospitals.

Specifications :

Technical Specifications

- **Chair Type:** 2-section dental chair with articulated headrest
 - **Operation:** Electrically controlled with foot control system
 - **Height Range:** 400 – 820 mm | **Backrest Tilt:** Up to 70°
 - **Upholstery:** Seamless, soft, replaceable, easy-to-clean
 - **Cuspidor:** Porcelain, swivel type
 - **Dentist Unit:** 2 air turbine + 1 micromotor, 3-in-1 syringe, ultrasonic scaler (25–30 kHz) with LED
 - **Control System:** Pneumatic / foot-controlled with coolant system
 - **Assistant Unit:** Saliva ejector + HVE suction, 3-in-1 syringe, independent controls
 - **Suction System:** Self-contained aspirator with auto cut-off protection
 - **Air Compressor:** Oil-free, noiseless, ¾ HP, ~100L tank, 5.6 kg/cm² pressure
 - **Dental Stools:** Adjustable height (500–700 mm), aluminium base with castors, lumbar support
 - **Digital X-Ray:** Portable handheld, 60 kV, 1–5 mA, 0.05–1.35 sec exposure, 16 GB storage, WiFi/Bluetooth, rechargeable battery
 - **Power Supply:** 240V, 50Hz
 - **Standard Accessories:** Handpieces set, ultrasonic scaler, intraoral X-ray with sensor, autoclave (19L), intraoral camera, LED curing light, dental instrument kit
-

