

Labycare-Ultra Low Freezers vertical -80° (SKSSULF384C)



Ultra Low Freezers vertical -80° manufacturer by **Labycare** in India. It is suitable for electronic industries, university laboratories, military industries, pelagic fishery companies.

Description :

Labycare Ultra Low Freezers vertical -80° range the internal volumes vary from 445 to 830 liter of safe and reliable storage capacity in a compact and user friendly design. Standard our freezers are equipped with stainless steel interior, inner doors and shelves. The combination of low power consumption, lowest noise level and the standard use of ecological, green refrigerants makes it the most environmental friendly ULT freezer available in the market.

Specifications :

Model	SKSSULF384C	SK
External Size (W*D*H)mm	1066*742*820	785
Internal Size (W*D*H)mm	450*439*514	440
Capacity	100 L	398
Type	Vertical	
Temp.Range	-10~ -80°C	
Temp. Accuracy	0.1 °C	
Control System	Microprocessor control,	Mic

		Large Screen LED Display		
Audio and Visual Alarm		High and low temperature, System blocking; Power failure		
Refrigeration Type		Direct refrigeration		
Refrigerant		Mixture refrigerant		
Compressor		1 set		
Condenser&Evaporator		Made of copper		
Construction	Structure	Unibody design, two-time foaming		
		100~155 mm foaming insulation,		
	Internal	SS304		
	Material			
	External Material	Cold-rolled steel coated with anti		
Door		Lockable door		
Caster		4 casters		
Chamber with Inner Door		2		4
Temp. Test Port		1-25 mm		2, 1
Consumption		1000 w		980
Power Supply		AC 110V/220V±10%,50/60HZ		
Standard Accessory		Temperature printer, record 7 US days data, 72 hours battery backup for printer and controller display		
Optional Accessory	Back-up System	CO2 Back-up system		
	Storage Racks	3 racks:4*3,3 racks:3*3		12
	and Boxes	63 Boxes		300
Package Size (W*D*H)mm		1210*870*1020		380
Gross Weight		140 kg		287

labycare

www.labycare.com | Email: labycare@gmail.com