

Labycare-Ultra Low Freezers Horizontal -80° (SKSS372C)



Ultra Low Horizontal -80° Deep Freezers are manufacturer by **Labycare** in India. These products are offered with thermostat or digital controller for effective monitoring of temperatures of these products. These products are manufactured using optimum quality raw materials that are sourced from reliable vendors. The inner walls of these products are insulated to avoid temperature loss.

Description :

Labycare Ultra Low Horizontal -80° Deep Freezers are easy to use and reduce cooling loss, saving energy. High efficiency air-cooled condenser, good heat dissipation and fast cooling. Ultra-thick environmental protection insulation layer could effectively reducing cooling loss, saving energy and prolonging service life. The inner body is made of stainless steel and the tube-on-sheet heat exchanger. Easy access to items, good cooling effect, easy to clean the liner. The control and the temperature indicating device are installed on the front of the box. You can adjust and check the temperature inside the box without leaving the door, beautiful, practical and generous. Microcomputer control, more precise and safe.

Specifications :

Model NO.

SKSSULF384D

Temperature Type

Single-temperature

Defrost Type

Manual Defrost

Function

Freezing

Certification

CE, RoHS

Name of Commodity

Freezer

Model

WI-86/120-Wa

Temp. Range

-40~-86 Degrees Celsius

Inner Capacity

120L

Power

500W

Shelves

0PCS

Inner Material

Stainless Steel

External Material

Steel Plates with Spraying

Compressor

Danfoss

Sensor

Ntc/PTC

Cooling Type

Air Cooling

Refrigerant

R404A

Power Supply

220V/50Hz; 110V/60Hz

Inner Dimension

500X400X560mm

Overall Dimension

1270X780X965mm

Trademark

WollenLab

Transport Package

Plywood Box

Specification

1270x780x965mm

Origin

INDIA

labycare

www.labycare.com | Email: labycare@gmail.com