

## Labycare-Ultra Chamber Refrigerator - 20° C (SKSSCR384F)



Ultra Chamber Refrigerator - 20° C manufacturer by **Labycare** in India.

### Description :

**Labycare Ultra Chamber Refrigerator - 20° C** is designed with two compressors and two chambers, which allows you to control the upper chamber and lower chamber independently. And the high-precision microcomputer temperature control system has the ability to set the temperature inside the cabinet in the range of -20°C ~ -40°C.

### Specifications :

- Capacity: 276 litres
- Ext. Dimensions (WxDxH mm): 920 x 805 x 1660
- Int. Dimensions (WxDxH mm): 600 x 480 x 970
- Temperature Range (°C): -20°C ~ -90°C
- **Insulation: 130 mm Non CFC, foamed in place polyurethane with VIP**
- Electrical: 110V, 220V 60Hz, 230V 50Hz
- Weight (kg): 256
- Refrigeration: 2 x 1.25 HP

---

labycare

www.labycare.com | Email: labycare@gmail.com