

## Labycare-Ultra Chamber Refrigerator Vertical -20°C (SKSSUCRV384G)



Ultra Chamber Refrigerator Vertical -20°C manufacturer by **Labycare** in India.

### **Description :**

**Labycare Low Temperature Freezer** offers a wide variety of research and storage application, such as low temperature scientific experiments, preservation of plasma, biomaterial, vaccine, biomedical properties of military products. It is suitable for research institutes, electronic industries, chemical engineering industries, hospitals, sanitation and antiepidemic stations, university laboratories, military industries.

### **Specifications :**

- External Size (W\*D\*H) mm: 675\*675\*1580
- Capacity: 268 L
- Type: Vertical
- Temperature Range: -20°C ~ -25°C
- Temperature: 0.1°C
- Control System: Microprocessor Control & LED Display
- Alarm: High and low temperature, System & sensor failure, Door ajar, Power failure alarm
- Refrigeration Type: Direct
- Refrigerant: R290
- Compressor: International famous compressor
- Condenser & Evaporator: High efficiency air-cooled condenser, made of copper
- Consumption: 593W
- Power Supply: AC 220V +-10%, 50/60Hz
- Gross Weight: 125 kg

---

labycare

www.labycare.com | Email: labycare@gmail.com