

## Labycare-Ultra Chamber Refrigerator Horizontal-20 (SKSSUCRH384H)



The horizontal system is the sample port on the back makes temperature monitoring convenient. The screen is easy to observe and writes a short record of the board. This will be passed to the next user green handle, making it easy to use 2 locks. Automatic cleaning of filter alarms and sensor error alarms. Adjustable shelf height, optional temperature recorder storage shelf and storage compartment.

### Description :

- "labycare" brings the ultra-high-tech Microprocessor controlled, solid-state Low Temperature Cabinet (Horizontal Models- Chest Opening type) micro control and CFC free refrigeration system. to meet various standard applications in diverse industrial applications such as medical, agricultural institutions, bio-chemical industrial and research laboratories.
- It is ideal laboratory freezer to provide a perfect freezing environment for the preservation of vaccines, blood plasma, test samples, specimens at desired temperature (Available in 3 Models i.e. -20°C, -30°C, -40°C.)

<b>Specifications :</b>	Horizontal Model No.	SKSSUCRH384H	SKSSUCRH384H	SKSSUCRH384H
Approx Volume (Ltr)		112	170	285
Capacity Cu.ft.		4 Cu.ft.	6 Cu.ft.	10 Cu.ft.
No. of Shelves		2	2	2
Chamber Size in MM		610 x 450 x 410	845 x 500 x 400	870 x 570 x 560
Temperature		Available in 3 Models i.e. -20°C, -30°C, -40°C.		
Temp. Controller		By Microprocessor Based PID Digital Temperature Indicator cum Controller.		
Temp. Accuracy		+1°C		
Display:		LED Display for Set Value(SV) and Process Value (PV)		
Relay		Solid state electronic relay with protective heat sink.		
Insulation		By High density PUF insulation		
Operations		Nearly silent operation with ultra-low vibration		
High Voltage Protection		Through automatic voltage stabilizer (Optional)		
Electric Supply		220/230V AC, 50/60Hz		

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)