

Labycare-Leak Test Apparatus. (SKSSLTA874)



Leak test apparatus This apparatus is an ideal instrument to check sealing integrity of Blister Packs; Strip Packs; Sachet and Pouch Packs etc. It is also used for checking leaks from Capped Bottles; Induction Sealed Containers; Sealed Or Crimped Tubes etc. The instrument is equipped with long lasting vacuum pump, which gives high level of vacuum (up to 550 Hg) in shortest possible time. Die pressed poly carbonate vacuum desiccator sustains vacuum for long time

Description :

Leak Tester works on the principle of creating a negative pressure in a closed chamber which inflates the sealed pouches due to pressure and leaked pouches don't inflate. To test the sample we place tablets sachet into the desiccator filled with Ethylene blue than vacuum the desiccators for a limited time. The vacuum gauge of the Leak test apparatus is connected to a control valve that isolates the vacuum inside desiccator from the vacuum pump.

Specifications :

labycare

www.labycare.com | Email: labycare@gmail.com