

Labycare-COD Digester. (SKSSCD876)



We are one of the most reliable companies in this domain and are into offering **COD Digester**.

Description :

COD Digesters is suitable for COD method for monitoring waste water in fluent and effluent, industrial process water and more. The instrument is used by Pollution Control Labs. Outer cabinet made of mild steel duly coated. Double wall construction with insulation for minimum heat losses increases the efficiency, of the system. Aluminum blocks with holes (Dimension 40 mm diameter X 80 mm depth) are used to hold 15 nos. COD tubes, Aluminum block are heated to maintain temp up to 180°C. Solid state digital temp controller, appropriate heater & insulation selection assure uniformity in all the samples. Selectable timer up to 120 mm with alarm is provided to set the digestion period (Refluxing). After digestion, analyze the sample with either colorimetric or titrimetric method.

Specifications :

Hole Size :

labycare

www.labycare.com | Email: labycare@gmail.com