

Labycare-Muffle Furnace -1200 (SKSSMF388C)



Muffle Furnace manufacturer by **Labycare** in India. It is box type heat treatment equipment used to change physical properties of samples at very high temperature; for example 1200°C, .

Description : **Labycare Muffle Furnace** is a laboratory furnaces are widely used in scientific experiments in physics lab, rice laboratories, steel and paint industries, biotech companies and small industrial production etc. Their major applications include general laboratory testing.

Specifications :	MAX WORKING	1200°C	VOLUME:	1 ~ 1000L
	TEMPERATURE:			
	HEATING ELEMENT:	High Quality Fe-Cr-Al-Mo Resistance Wire	FURNACE DOOR:	The Upper Door Can Be Customized Away From The Main Body With Over Temp Protection
	TEMPERATURE CONTROL:	50 Segments Programmable And Auto Control	PIDALARM:	
	COOLING:	Cooling Fan Can Reduce The Temperature Of The Casing	VOTAGE:	AC 110V Or AC 220V
	MAX. HEATING RATE:	20°C /min	TEMPERATURE ACCURACY:	±1°
			APPLICATION:	Industrial Temperature Treatment