

Labycare-Horizontal Laminar Air Flow (SKSSHLAF394C)



Horizontal Laminar Air flow manufacturer by **Labycare** in India. It provided by us facilitates easier installation and has lower needs for maintenance. This bench is acclaimed for its sturdiness.

Description :

The cabinet are fabricated out of thick mild steel cabinet duly powdered coated. The work table is made of thick stainless steel sheet. Side panels are provided of thick transparent acrylic sheets. The unit is fitted with pre-filter and is made to pass through highly efficiency 99.9%, thus retaining all air borne particles of size: 0.3 micron and larger. Using a dynamic machine, the blower and motor assembly is statically and dynamically balanced. Blowers are mounted on vibration isolation mounts to minimize vibration. The working is illuminated by fluorescent lighting fitted to the unit. Height of the working table provides comfortable "Sit-Down" working position for the operating to work on 220/230 volts AC supply.

Specifications :

- Mode of 19 mm ply-board & stainless steel.
- Epoxy painted inside.
- HEPA filter efficiency 99.97% upto 0.3 micron.
- Protection grill for HEPA filter.
- Pre-filter efficiency 99% upto 5 microns.
- Stainless Steel work bench.
- UV germicidal lamp.
- Fluorescent illumination.
- Cock for gas.
- Static pressure Manometer.
- Balanced blower with motor.

labycare

www.labycare.com | Email: labycare@gmail.com