

Labycare-Pass Box Dynamic (SKSS399)



Pass Box Dynamic manufacturer by Labycare in India.

Description :

Labycare Pass Box Dynamic is a kind of class 100 purification equipment for Pharmaceutical clean room. One hand, it plays the role of the air lock, Does not let high grade rank clean room pressure be leaked. On the other hand, the self purification effect is realized in the process of pass box to ensure the clean of the articles entering the clean area and reduce the cross pollution caused by the objects. Customized design is available. Pass box inside bottom corner is round shape so that there is any dust accumulation External size: 800*670*1150 mm Internal size: 600*600*600 mm External Material: stainless steel304 (t=1.2) Inside Material: stainless steel304 (t=1.2) Double layer vacuum glass Both side door is inlaid type Locker: software computer interlocking DOP Gel Sealant Mini-Pleats HEPA air filter Air speed: 0.45 m/s adjustable DOP/POA port: 1set Introduce: Dynamic Pass Box.

Specifications :

- Interlocking: Mechanical interlocking for the doors.
- UV Light: UV light with hour meter with interlocking arrangement to put-off UV.
- Construction: Floor mounted construction with three side wall structure & door drop seal.
- Support Stand: Support stand to mount the pass box.
- Flange: Flange for sealing the gap between the pass box and the clean room wall.
- Door: Swing type doors with view panels.
- Motor-Blower: Statically & Dynamically balanced motor-blower with suspension arrangement to reduce noise level.
- Door open indicator.
- Optional UV Lamp & Fluorescent tube light.

labycare

www.labycare.com | Email: labycare@gmail.com