

## Labycare-Magnetic Stirrer with Hot Plate (SKSS346)



This group of magnetic stirrers is similar to the conventional magnetic stirrers but have additional hot plate to achieve heating of solutions to be stirred. The heating is controlled by energy regulator, which ensures regulation at any desired percentage of the capacity of heater.

### Description :

**SKSS** introduced world's best Magnetic-stirrers. This group of magnetic stirrers is similar to the conventional magnetic stirrers but have additional hot plate to achieve heating of solutions to be stirred. The heating is controlled by energy regulator, which ensures regulation at any desired percentage of the capacity of heater.

### Specifications :

- .Maximum Speed RPM: 1800
- .Stirring Capacity Ltrs: 2liter
- .Heating Capacity Watts 500
- .Dimension ( W x D x H ) mm 150 x150 x 170
- Stirring Paddle ØxL mm: 2

#### • FEATURES

- Aluminum with ceramic coating hotplate, anti-corrosion and durable.
- With strong magnetic to prevent the magnetic stir bar from escaping.
- Stepless speed regulation, running smoothly at low speed and efficiently at high speed.
- Over temperature protection with independent safety circuit design, to ensure the lab safety.

- Knob design for easy to control temperature and speed (ATO-HS-19 only).
  - Digital hotplate magnetic stirrer with LCD display, easy to read (ATO-HSC-19T only).
  - Adjustable sensor fixing bracket, convenient to use.
- 

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)