

Labycare-Laboratory Refrigerator (SKSSLR369B)



Laboratory Refrigerators manufacturer by **Labycare** in India. It supplied with digital controller, glass door, alarm system to meet the demanding requirements of laboratory research. These refrigerator has temperature ranging from +2°C to +8°C and all models are combined with two exact sensors and auto defrost.

Description :

- **Labycare laboratory fridge** is used by healthcare facilities and professionals to store samples and specimens, as well as vaccines and medicines, at a specific temperature to ensure they do not become spoiled.

Specifications :

• Model Number	SKLR369B
• Capacity	360L
• Display	Digital
• Brand	Labycare
• Net/Gross(kg)	103/112 kg
• Exterior Dimension(mm)(WxDxH)	620 x 655 x 1955mm
• Interior Dimension(mm)(WxDxH)	555 x 460 x 1340mm
• Temprange	2 to 8 deg C
• Voltage	220V
• Input power(w)	430W
• Frequency	50 Hz

labycare

www.labycare.com | Email: labycare@gmail.com