

Labycare-Elisa Washer. (SKSS394)



ELISA test kits with micro-titre plates, ELISA washer , plate incubator
UPS of sufficient capacity to ensure uninterrupted finishing of ongoing testing

Description :

A washer for microtitre plates, designed to ensure thorough washing of reagents between ELISA steps. The equipment is used for many types of ELISA.

To be used by appropriately trained staff.

Often used in reference or referral testing and in high throughput screening.

Manufacturers recommended waste disposal:

Specifications :

Micro-plate or strip washer for ELISA (Enzyme-Linked Immunosorbent Assay)

Platform: bench-top, open and automatic system

8-channel strip manifold

Accommodates flat, round, U or V-shaped bottom plates

75 user-definable protocols

Wash parameters include: 16-character assay name, number of cycles, wash volume, flow rate and variable soak times, aspiration and dispensing

Crosswise aspiration/double aspiration of flat bottom micro-plates for reduced residual liquid

Built-in vacuum & pressure pump assembly

With bottles for waste, rinse and wash solutions

Wash cycles: Between 1 to 10

Dispensing volumes from 25 to 3000ul

Soak time 1-600 seconds

Fluid flow rate in 150 to 1000ul/well/second to accommodate cellular assays

Spill-over protection & electronics isolated from fluidics

labycare

www.labycare.com | Email: labycare@gmail.com