

## Labycare-RF Blood Tube Sealer (SKSS388)



**Blood Bank Tube Sealer** manufacturer by **Labycare** in India. It is required to seal the blood using high frequency dielectric heating mechanism. It is the safest way to preserve the blood while it can be examined later on demand. The process requires that the machine must be free from any defects. It is also necessary that the equipment operates under the preferred conditions to make sure that there is no mistake and the blood is not contaminated.

### Description :

With thorough understanding of the relevant industry, we are engaged in providing a unique range of **Labycare Blood Bag Tube Sealer**. Offered tube sealer is designed to seal tubing of PVC bags and provide hermetic seal on plastic tubing. This tube sealer is developed by quality tested components and contemporary technology by our adroit professionals. Our provided tube sealer has got high appreciation owing to its automatically triggers sealing.

Parameter	WeldOn – B
Dimension (W x D x H) (mm)	220 X 200 X 120
Weight (Kg)	4
RF Frequency	40.68 MHZ
RF Output	80 watt
Seal time	1.5 Sec (auto adjustable)
Sealing Type	Automatic
EEC Standard	Class 1 (93/42 and 47/2007 MDD)
Input Voltage	115/230 V, 50/60 Hz
Ref Standard	E 60601-1 for electric safety, E 60601-1-2 for EMC

labycare

www.labycare.com | Email: labycare@gmail.com